

Laboratory Freezers



Sparkfree Laboratory Freezers
Advanced Laboratory Freezers
Sparkfree Laboratory Fridge/Freezers

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Sparkfree Laboratory Freezers

The Labcold Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Everything that could conceivably cause a spark, such as the controller, is located outside of the freezer chamber making these freezers an important part of laboratory safety. Designed specifically for laboratories, these freezers are solidly built and come in both upright and chest configuration. All Labcold Sparkfree freezers and chest freezers are factory fitted with a door lock for optimum security.



[Labcold RLVF02203 Sparkfree Laboratory Freezer 36 Litres](#)

Capacity: 36 Litre
Config: Upright Freezer
515mm H x 440mm W x 500mm D
Product Code: FPD-09917



[Labcold RLVF02203/LK Sparkfree Laboratory Freezer 36 Litres Lockable](#)

Capacity: 36 Litre
Config: Upright Freezer
515mm H x 440mm W x 500mm D
Product Code: FPD-09918



[LEC ISU27C Laboratory Freezer 50 Litres](#)

Capacity: 50 Litre
Config: Bench Top
562mm H x 500mm W x 510mm D
Product Code: FPD-09804



[LEC ISU37C Laboratory Freezer 102 Litres](#)

Capacity: 102 Litre
Config: Under Bench
830mm H x 495mm W x 620mm D
Product Code: FPD-09805



[Labcold RLVF04202A Sparkfree Laboratory Freezer 120 Litres](#)

Capacity: 120 Litre
Config: Upright Freezer
850mm H x 545mm W x 565mm D
Product Code: FPD-07539



[Labcold RLVF04202A/LK Sparkfree Laboratory Freezer 120 Litres](#)

Capacity: 120 Litre
Config: Upright Freezer
850mm H x 545mm W x 565mm D
Product Code: FPD-07540



[Liebherr LGUex 1500 Laboratory Freezer - 139 Litre](#)

Capacity: 139 Litre
Config: Bench Top
820mm H x 600mm W x 615mm D
Product Code: FPD-10261



[LEC LSF151UK Laboratory Freezer 151 Litres](#)

Capacity: 151 Litre
Config: Under Bench
845mm H x 595mm W x 595mm D
Product Code: FPD-09806

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Sparkfree Laboratory Freezers

The Labcold Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Everything that could conceivably cause a spark, such as the controller, is located outside of the freezer chamber making these freezers an important part of laboratory safety. Designed specifically for laboratories, these freezers are solidly built and come in both upright and chest configuration. All Labcold Sparkfree freezers and chest freezers are factory fitted with a door lock for optimum security.



[Labcold RLVF07203 - 180 Litre Sparkfree Laboratory Freezer with Solid Door](#)

Capacity: 180 Litre
Config: Upright Freezer
1440mm H x 544mm W x 570mm D
Product Code: FPD-07708



[Labcold RLVF07203LK - 180 Litre Sparkfree Laboratory Freezer with Solid Door](#)

Capacity: 180 Litre
Config: Upright Freezer
1440mm H x 544mm W x 570mm D
Product Code: FPD-07709



[Labcold RLCF0720 Sparkfree Laboratory Freezer 215 Litres](#)

Capacity: 215 Litre
Config: Chest Freezer
860mm H x 720mm W x 730mm D
Product Code: FPD-04934



[LEC LSF232UK Laboratory Freezer 232 Litres](#)

Capacity: 232 Litre
Config: Free Standing
1570mm H x 595mm W x 634mm D
Product Code: FPD-09807



[Labcold RLVF14201 Sparkfree Laboratory Freezer 375 Litres](#)

Capacity: 375 Litre
Config: Upright Freezer
1850mm H x 595mm W x 640mm D
Product Code: FPD-04933



[Labcold RLCF1520 Sparkfree Laboratory Freezer 447 Litres](#)

Capacity: 447 Litre
Config: Chest Freezer
860mm H x 1300mm W x 730mm D
Product Code: FPD-04935



[Labcold RLCF2120 Sparkfree Laboratory Freezer 607 Litres](#)

Capacity: 607 Litre
Config: Chest Freezer
860mm H x 1700mm W x 730mm D
Product Code: FPD-04936

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF02203 Sparkfree Laboratory Freezer 36 Litres

Product Code: FPD-09917

Dimensions: 515mm H x 440mm W x 500mm D



- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future
- Capacity: 36 Litres (gross)
- Construction: White enamel exterior, fully moulded plastic interior
- Shelves supplied: 1, wire mesh
- Sparkfree interior
- Counter- top use
- CFC and HCFC free
- Temperature set-point range: -10 to -26oC
- Operator controlled defrost

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC ISU27C Laboratory Freezer 50 Litres



Product Code: FPD-09804

Dimensions: 562mm H x 500mm W x 510mm D

Internal Dimensions: 365mm H x 345mm W x 330mm D

- -18°C to -20°C set temperature - Perfect for short term sample storage when space is at a premium
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark-free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- Small countertop size - For limited space needing biomedical and laboratory refrigeration
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Colour: White/green
- Door Type: Solid
- Door Style: Hinged, reversible
- Door Open Depth: 976mm
- Lock: Yes
- One fixed wire shelf
- Manual defrost
- Climate Class: N
- Refrigerant Type: R600a
- CFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Noise Level: <45dB
- Insulation Thickness: 30mm
- Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC ISU37C Laboratory Freezer 102 Litres



Product Code: FPD-09805

Dimensions: 830mm H x 495mm W x 620mm D

Internal Dimensions: 683mm H x 414mm W x 415mm D

- -18°C to -20°C set temperature - Perfect for short term sample storage when space is at a premium
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Colour: White/green
- Door Type: Solid, Hinged-reversible
- Door Open Depth: 1050mm
- Lock: Yes
- Adjustable Feet: 1, front
- Two adjustable wire shelves
- Manual defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free & HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Noise Level: <43dB
- Insulation Thickness: 35mm
- Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 1.3

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF04202A/LK Sparkfree Laboratory Freezer 120 Litres



Product Code: FPD-07540

Dimensions: 850mm H x 545mm W x 565mm D

- Capacity 120 Litres
- Configuration: Under-bench lockable
- Sparkfree safety interior
- CFC, HCFC and Amonia free
- Operating temperature -20°C
- Operator controlled defrost
- Maximum ambient temperatures up to 32°C
- Low energy
- Door lock
- 4 Shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temperature display
- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF04202A Sparkfree Laboratory Freezer 120 Litres



Product Code: FPD-07539

Dimensions: 850mm H x 545mm W x 565mm D

- Capacity 120 Litres
- Configuration: Under-bench
- Sparkfree safety interior
- CFC, HCFC and Amonia free
- Operating temperature -20°C
- Operator controlled defrost
- Maximum ambient temperatures up to 32°C
- Low energy
- 4 Plastic removeable drawers separated by fixed refrigerated shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temperature display
- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Liebherr LGUex 1500 Laboratory Freezer - 139 Litre



Product Code: FPD-10261

Dimensions: 820mm H x 600mm W x 615mm D

- Temperature range -9°C to -26°C
- Suitable for the storage of explosive or highly inflammable substances
- Digital Temperature Control & Display
- Temperature deviation alarm with volt free contacts
- Door open alarm
- Manual Defrost
- Door Lock and 2 keys
- CFC & Ammonia free
- 3 x drawers + 1 Basket
- 220-230v 50Hz
- Max. Ambient +38°C

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP_{TM}

LEC LSF151UK Laboratory Freezer 151 Litres



Product Code: FPD-09806

Dimensions: 845mm H x 595mm W x 595mm D

- -18°C to -25°C set temperature - Perfect for short term sample storage when space it at a premium
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Roller castors - To make the freezer easier to move for cleaning and maintenance
- External Finish: White, Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Open Depth: 1165mm
- Open Door Sensor
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Three fixed wire basket Shelves: 3
- Manual defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Aux Output
- Insulation Thickness: 45mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF07203LK - 180 Litre Sparkfree Laboratory Freezer with Solid Door



Product Code: FPD-07709

Dimensions: 1440mm H x 544mm W x 570mm D

- Capacity: 180 Litres (gross)
- Construction: White enamel exterior, plastic interior
- Door Type: Solid
- Internal Arrangement: 4 drawers and 2 compartments with pull down flaps separated by fixed refrigerated shelves
- Sparkfree Interior
- CFC and HCFC free
- Temperatures set-point range: -18oC to -26oC
- Analogue thermostat control
- Visual alarm for high temperature
- Low energy operation
- Operator controlled defrost
- Locking Door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF07203 - 180 Litre Sparkfree Laboratory Freezer with Solid Door



Product Code: FPD-07708

Dimensions: 1440mm H x 544mm W x 570mm D

- Capacity: 180 Litres (gross)
- Construction: White enamel exterior, plastic interior
- Door Type: Solid
- Internal Arrangement: 4 drawers and 2 compartments with pull down flaps separated by fixed refrigerated shelves
- Sparkfree Interior
- CFC and HCFC free
- Temperatures set-point range: -18oC to -26oC
- Analogue thermostat control
- Visual alarm for high temperature
- Low energy operation
- Operator controlled defrost
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLCF0720 Sparkfree Laboratory Freezer 215 Litres



Product Code: FPD-04934

Dimensions: 860mm H x 720mm W x 730mm D

- Configuration: Chest
- Capacity 215 Litres
- Operating temperature -20°C
- Sparkfree interior
- Door lock
- Alarm
- Low energy
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LSF232UK Laboratory Freezer 232 Litres



Product Code: FPD-09807

Dimensions: 1570mm H x 595mm W x 634mm D

- -18°C to -25°C Set temperature - Perfect for short term storage for a large amount of samples
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark-free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Roller castors - To make the freezer easier to move for cleaning and maintenance
- External Finish: White, Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Open Depth: 1185mm
- Open Door Sensor
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Four wire drawers, 2 fixed flaps, Four Moulded Drawer
- Manual defrost
- Climate Class: SN, Refrigerant Type: R600a
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Aux Output: Yes
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF14201 Sparkfree Laboratory Freezer 375 Litres



Product Code: FPD-04933

Dimensions: 1850mm H x 595mm W x 640mm D

- Capacity 375 Litres
- Configuration: Free-standing
- Sparkfree safety interior
- CFC, HCFC and Amonia free
- Operating temperature -20°C
- Operator controlled defrost
- Maximum ambient temperatures up to 32°C
- 6 Shelves
- Temperature Display
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLCF1520 Sparkfree Laboratory Freezer 447 Litres



Product Code: FPD-04935

Dimensions: 860mm H x 1300mm W x 730mm D

- Configuration: Chest
- Capacity 447 Litres
- Operating temperature -20°C
- Sparkfree interior
- Door lock
- Alarm
- Low energy
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLCF2120 Sparkfree Laboratory Freezer 607 Litres



Product Code: FPD-04936

Dimensions: 860mm H x 1700mm W x 730mm D

- Configuration: Chest
- Capacity 607 Litres
- Operating temperature -20°C
- Sparkfree interior
- Door lock
- Alarm
- Low energy
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Advanced Laboratory Freezers

Advanced freezers are completely constructed from powder coated mild steel with a stainless steel interior so are built to last. Insulated with Envirofoam and fitted with micro processor controllers these freezers provide accurate temperature control with low GWP. All Advanced freezers have optimised air circulation for superior temperature control and feature a temperature display, high and low temperature alarms and door locks fitted as standard. For ease of positioning the larger cabinets have heavy duty wheels with a brake. Designed for any situation where accurate temperature control is vital and there is a need to store a large quantity of product. Advanced freezers can be fitted with access ports and remote alarm contacts.



[Labcold RAFR05202 Advanced Laboratory Freezer 150 Litres](#)

Capacity: 150 Litre
Config: Upright Freezer
830mm H x 600mm W x 620mm D
Product Code: FPD-04946



[Labcold RAFR21262 Advanced Laboratory Freezer 571 Litres](#)

Capacity: 571 Litre
Config: Upright Freezer
1980mm H x 720mm W x 820mm D
Product Code: FPD-04948



[Labcold RAFR44262 Advanced Laboratory Freezer 1245 Litres](#)

Capacity: 1245 Litre
Config: Upright Freezer
1980mm H x 1440mm W x 820mm D
Product Code: FPD-04949

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RFR05202 Advanced Laboratory Freezer 150 Litres



Product Code: FPD-04946

Dimensions: 830mm H x 600mm W x 620mm D

- Configuration: Underbench, solid door
- Capacity: 150 Litres
- Operating temperature: -20°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 3 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RAFR44262 Advanced Laboratory Freezer 1245 Litres



Product Code: FPD-04949

Dimensions: 1980mm H x 1440mm W x 820mm D

- Configuration: Free-standing, double solid door
- Capacity: 1245 Litres
- Operating temperature: -20°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 8 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RFR21262 Advanced Laboratory Freezer 571 Litres



Product Code: FPD-04948

Dimensions: 1980mm H x 720mm W x 820mm D

- Configuration: Free-standing, solid door
- Capacity: 571 Litres
- Operating temperature: -20°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 4 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Combined Laboratory Fridge/Freezers

This is the latest generation of combined sparkfree laboratory refrigerator and freezer. Providing all the routine storage required on a floor area no bigger than a standard fridge. For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future. Both refrigerated and frozen samples can be stored and accessed within the same small floor area.



[LEC ISR27C Laboratory Fridge Freezer](#)

Capacity: 83 Litre
Config: Bench Top
660mm H x 502mm W x 540mm D
Product Code: FPD-09808



[LEC LSC119UK Laboratory Fridge Freezer](#)

Capacity: 119 Litre
Config: Under Bench
860mm H x 540mm W x 600mm D
Product Code: FPD-09809



[LEC LSC263UK Laboratory Fridge Freezer](#)

Capacity: 263 Litre
Config: Upright Refrigerator
1680mm H x 540mm W x 600mm D
Product Code: FPD-09810



[LEC LSC324UK Laboratory Fridge Freezer](#)

Capacity: 324 Litre
Config: Upright Refrigerator
1870mm H x 600mm W x 600mm D
Product Code: FPD-09811



[Labcold RLFF13246 Laboratory Fridge/Freezer 230/160 Litres](#)

Capacity: 390 Litre
1965mm H x 595mm W x 630mm D
Product Code: FPD-04937



[Labcold RLFF13246/LK Laboratory Fridge/Freezer 230/160 Litres](#)

Capacity: 390 Litre
1965mm H x 595mm W x 630mm D
Product Code: FPD-04938

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF02203/LK Sparkfree Laboratory Freezer 36 Litres Lockable

Product Code: FPD-09918

Dimensions: 515mm H x 440mm W x 500mm D



- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future
- Capacity: 36 Litres (gross)
- Construction: White enamel exterior, fully moulded plastic interior
- Shelves supplied: 1, wire mesh
- Sparkfree interior
- Counter- top use
- CFC and HCFC free
- Temperature set-point range: -10 to -26oC
- Operator controlled defrost

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC ISR27C Laboratory Fridge Freezer



Product Code: FPD-09808

Dimensions: 660mm H x 502mm W x 540mm D

- +2°C to +10°C fridge set temperature, -6°C to -9°C freezer set temperature - Designed for when both refrigeration and freezing is needed, perfect for smaller spaces
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge and freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the fridge or freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- Small countertop size - For limited space needing biomedical and laboratory refrigeration
- Fridge Net Capacity: 75, Freezer Net Capacity: 5.5
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Colour: White/green
- Door Type: Solid
- Door Open Depth: 976mm
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Two wire Shelves, Shelf Size: 438x274x6mm & 438x137x6mm
- Fridge: Automatic defrost, Freezer: Manual defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, CFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Noise Level: <43dB
- Insulation Thickness: 30mm, Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 1.74

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LSC119UK Laboratory Fridge Freezer



Product Code: FPD-09809

Dimensions: 860mm H x 540mm W x 600mm D

- Set temperature - Designed for when both refrigeration and freezing is needed, perfect for smaller spaces
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge and freezer
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the fridge-freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Roller castors - To make the fridge-freezer easier to move for cleaning and maintenance
- Fridge Net Capacity: 105L, Freezer Net Capacity: 12L
- External Finish: White, Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Temperature Range: Fridge +2°C to +10°C, freezer -18°C to -25°C
- Colour: White/green
- Door Type: Solid, Door Open Depth: 1090mm, Open Door Sensor
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Two wire shelves & One solid moulded Drawer, Two full adjustable door bins + two half
- Fridge: Automatic defrost, Freezer: Manual defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Aux Output: Yes
- Insulation Thickness: 45mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LSC263UK Laboratory Fridge Freezer



Product Code: FPD-09810

Dimensions: 1680mm H x 540mm W x 600mm D

- Our spark free laboratory fridge-freezer offers great flexibility and is ideal for the safe storage of a wide variety of volatile liquids and chemicals.
- Set temperature - Designed for when both refrigeration and freezing is needed, perfect for smaller spaces
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge and freezer
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the fridge-freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Reversible door - For maximum versatility in sometimes limited working spaces
- Roller castors - To make the fridge-freezer easier to move for cleaning and maintenance
- External Finish: White, Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Temperature Range: Fridge +2°C to +10°C, freezer -18°C to -25°C
- Door Type: Solid, Door Open Depth: 1110mm, Open Door Sensor
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Three wire Shelves: 3 (fridge), Moulded Drawer: Fridge: 2 solid Adjustable Door Bins: 2 full + 2 half
- Fridge: Automatic defrost, Freezer: Manual defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane
- Fridge Net Capacity: 157 Litre, Freezer Net Capacity: 88 Litre

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LSC324UK Laboratory Fridge Freezer



Product Code: FPD-09811

Dimensions: 1870mm H x 600mm W x 600mm D

- Our new LSC324UK is our largest fridge-freezer unit yet, offering you even greater storage capacity
- Set temperature - Designed for when both refrigeration and freezing is needed, perfect for smaller spaces
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge and freezer
- Antimicrobial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark-free - There is nothing electrical inside the fridge-freezer, so it is safe to store reactive substances in without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Reversible door - for maximum versatility in sometimes limited working spaces
- Roller castors - To make the fridge-freezer easier to move for cleaning and maintenance
- External Finish: White, Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Temperature Range: Fridge +2°C to +10°C, freezer -18°C to -25°C
- Door Type: Solid, Door Open Depth: 1190mm, Open Door Sensor
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Three wire Shelves, Two solid moulded drawers, Adjustable Door Bins: 2 full + 4 half
- Fridge: Automatic defrost, Freezer: Manual defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLFF13246/LK Laboratory Fridge/Freezer 230/160 Litres



Product Code: FPD-04938

Dimensions: 1965mm H x 595mm W x 630mm D

- Lab Fridge Freezer
- Configuration: Freestanding
- Capacity 170/160 Litres
- +4/-20°C
- Door lock
- Sparkfree interior
- Auto defrosting on freezer
- Low energy
- 4 shelves in refrigeration
- 5 shelves in freezer
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLFF13246 Laboratory Fridge/Freezer 230/160 Litres



Product Code: FPD-04937

Dimensions: 1965mm H x 595mm W x 630mm D

- Lab Fridge Freezer
- Configuration: Freestanding
- Capacity 170/160 Litres
- +4/-20°C
- Sparkfree interior
- Auto defrosting on freezer
- Low energy
- 4 shelves in refrigeration
- 5 shelves in freezer
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696

Fax 0115 727 0123

Email info@pharmacy-fridge.co.uk

How To Order

NHS Purchase Order, MOD Purchase Order, Local Authority Purchase Order

All official government, NHS, MOD, Police, Educational and Local Authority purchase orders are accepted. Fax your purchase order to 0115 727 0123 or email to info@pharmacy-fridge.co.uk.

Bank Transfer

Sent bank transfers to:

Payee	FPD Group Limited		
Bank	HSBC Bank plc		
Sortcode	401903	IBAN	GB97MIDL40190302074613
Account Number	02074613	Swift	MIDLGB2104A

If possible send confirmation of your transaction, to speed up the order process.

Cheques

Make cheques payable to FPD Group Limited and post to:

FPD Group Limited
The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Credit/Debit Card

Payment from UK based businesses can be taken over the telephone or processed via the website [shopping cart](#). Customers with issued quotations may process their order online using our [Purchase Code](#) web page.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Delivery

48 Hour delivery is available on most refrigerators through-out the UK mainland, larger units may require a longer lead time, please contact our sales staff whilst ordering.

Approximate Delivery lead times can be confirmed during our quotation process, typical delivery time is 10-14 days from receipt of payment.

**Drivers will not carry goods into buildings.

Goods normally arrive by truck on a pallet and are delivered Kerbside. Please note you may need lifting equipment or a number of strong workers to handle the order when it arrives. If you have any special requirements please [contact us](#).

Report any damages or shortages immediately, please read delivery notes carefully and use contact numbers provided or call FPD Group Ltd.

For more detail refer to [FPD Group Ltd Terms and Conditions](#) documentation.

**Large Products will require a two person delivery and installation, this is built into the price, the larger fridges also may require a longer delivery time of 3-4 weeks please check whilst ordering.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP[™]

Pharmacy Fridge Exports

Since 1961 FPD Group Ltd's family have been supplying medical equipment. Our extensive experience and technical knowledge now concentrates on providing clients with a quality selection of products for hospitals and healthcare organisations worldwide. We receive export enquiries and orders from all parts of the world and have established customers in UK, Europe, Africa and Asia. Export business activities are supervised by either our UK or Asia offices.

Our client base ranges from overseas ministry of health, defence, hospitals, surgeries, clinics and doctors to aid agencies and disaster relief. FPD Group Ltd export almost any medical product listed on our websites.

FPD estimators have the expertise to supply a single item or to fully equip many types of healthcare facilities. As well as the initial supply of equipment we also offer room design services using our in-house Autodesk Revit Architecture 2013 CAD systems.

** Export items are normally subject to minimum order values, additional protective packing charges and export shipping costs. Prices shown are Ex Works unless otherwise agreed in writing.*



FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

About Pharmacy Fridge

Pharmacy Fridge supply a wide range of fridges meeting the specifications of the healthcare industry. To buy simply add your required product to the shopping cart and submit your order.

For dispensary, surgery, ward and general use for the storage of drugs and vaccines under controlled conditions in line with the "Code of Ethics and Standards" of the Royal Pharmaceutical Society of Great Britain. Where necessary to achieve the required temperature control internal air-circulating fans are fitted.

Temperature display, alarms, max/min temperature recording, and door locks are standard and all models are CFC, HCFC and Ammonia free.

Are you storing your vaccines correctly? Recent research has shown an alarming level of use of domestic refrigerators to store vaccines. Did you know that not only can this damage the effectiveness of the vaccines it also poses a security risk for medicines that should be locked away?

FPD Group Ltd

FPD Group Limited is a family business with over 50 years experience in the supply and installation of a wide range of Healthcare Equipment, Medical Refrigeration and Hygienic Storage Furniture. Our customer base includes organisations such as NHS Trusts, PCT's, HMP Prison service, Ministry of Defence, GP Practices, Pharmacy Shops and Dispensaries, Universities, Schools, Dental Practices, Laboratories, Veterinary Practices, Care Homes and Private Buyers. We also export to the Irish Republic, Asia, United Arab Emirates and Qatar. In the UK the Company is based in Exeter, Devon with an additional office in Asia.

Our product range includes pharmacy storage solutions, refrigeration, controlled drugs storage, medical trolleys, clinical waste bins, medical lighting solutions, medical and healthcare furniture ranges and a number of bespoke solutions on request.

FPD Group has an in-house CAD visualisation department and accepts most file formats from your architects' plans, please ask for a quotation.

We are fully committed to the provision of a comprehensive range of professional products, services and furniture. Our philosophy is to provide a complete and appropriate solution to our customers' requirements.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Contact Pharmacy Fridge

FPD Group Ltd have over 40 years experience in the healthcare and shopfitting industry. We offer consultancy, design and installation services. Along with our Pharmacy Shelving and Storage solutions we also provide a wide range of medical equipment and healthcare furniture. Please visit www.fpd-group.co.uk for more information on our medical products.

For general enquiries complete the form below and we will contact you ASAP, please provide as much information as possible including dimensions where required. Please email any plans or photos of installations directly to sales@medicine-fridge.co.uk.

Name	<input type="text"/>	Address	<input type="text"/>
Company	<input type="text"/>	...	<input type="text"/>
Telephone	<input type="text"/>	City	<input type="text"/>
Fax	<input type="text"/>	State	<input type="text"/>
Email	<input type="text"/>	PostalCode	<input type="text"/>
		Country	<input type="text"/>
Comments	<input type="text"/>		

Submit Enquiry

Sales

Tel: 0113 350 8696
Fax: 0115 727 0123
Email: sales@pharmacy-shelving.co.uk

Payments & Billing

Tel: 01392 927 567

Address

FPD Group Limited
The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk