

Pharmacy & Vaccine Refrigeration



Pharmacy Refrigeration
LEC Medical
Labcold
Shoreline Medical
Dometic Silenco
Advanced Pharmacy Refrigeration

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[28 Litre Silenco Lockable Pharmacy Absorption Fridge](#)

£430.00

Capacity: 28 Litre
Config: Bench Top
580mm H x 422mm W x 393mm D

Product Code: FPD-09818



[Labcold IntelliCold™ RLDF0110 - 36 Litre Pharmacy Fridge with Solid Door](#)

£495.00

Capacity: 36 Litre
Config: Bench Top/Wall Mounted
535mm H x 445mm W x 465mm D

Product Code: FPD-09835



[LEC PG102C - Pharmacy Fridge 41 Litre with Glass Door](#)

£695.00

Capacity: 41 Litre
Config: Bench Top
702mm H x 379mm W x 500mm D

Product Code: FPD-09752



[LEC PE102C - Pharmacy Fridge 41 Litre with Solid Door](#)

£585.00

Capacity: 41 Litre
Config: Bench Top
702mm H x 379mm W x 500mm D

Product Code: FPD-09742



[LEC PE109C - Pharmacy Fridge 45 Litre with Solid Door](#)

£441.00

Capacity: 45 Litre
Config: Bench Top
530mm H x 470mm W x 465mm D

Product Code: FPD-09743



[Shoreline SM45 - 45 Litres Pharmacy Fridge with Solid Door](#)

£478.00

Capacity: 45 Litre
Config: Bench Top
715mm H x 380mm W x 375mm D

Product Code: FPD-06627



[53 Litre Silenco Lockable Pharmacy Absorption Fridge](#)

£515.00

Capacity: 53 Litre
Config: Bench Top
592mm H x 486mm W x 494mm D

Product Code: FPD-09819



[Shoreline SM60G - 60 Litres Pharmacy Fridge with Glass Door](#)

£718.00

Capacity: 60 Litre
Config: Bench Top
660mm H x 470mm W x 455mm D

Product Code: FPD-04358

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[Labcold IntelliCold™ RLDF0210 - 66 Litre Pharmacy Fridge with Solid Door](#)

£590.00

Capacity: 66 Litre
Config: Bench Top/Wall Mounted
725mm H x 450mm W x 485mm D

Product Code: FPD-09836



[Shoreline SM75 - 75 Litres Pharmacy Fridge with Solid Door](#)

£608.00

Capacity: 75 Litre
Config: Bench Top
790mm H x 450mm W x 430mm D

Product Code: FPD-07731



[LEC PE207C - Pharmacy Fridge 82 Litre with Solid Door](#)

£754.00

Capacity: 82 Litre
Config: Bench Top
660mm H x 502mm W x 540mm D

Product Code: FPD-09744



[LEC PG207C - Pharmacy Fridge 82 Litre with Glass Door](#)

£999.00

Capacity: 82 Litre
Config: Bench Top
660mm H x 502mm W x 450mm D

Product Code: FPD-09753



[LEC PG307C - Pharmacy Fridge 107 Litre with Glass Door](#)

£1,213.00

Capacity: 107 Litre
Config: Under Bench
830mm H x 495mm W x 620mm D

Product Code: FPD-09754



[LEC PE307C - Pharmacy Fridge 107 Litre with Solid Door](#)

£980.00

Capacity: 107 Litre
Config: Under Bench
830mm H x 495mm W x 620mm D

Product Code: FPD-09745



[Shoreline SM161G - 141 Litres Pharmacy Fridge with Glass Door](#)

£1,128.00

Capacity: 141 Litre
Config: Under Bench
835mm H x 600mm W x 615mm D

Product Code: FPD-08130



[Shoreline SM161 - 141 Litres Pharmacy Fridge with Solid Door](#)

£938.00

Capacity: 141 Litre
Config: Under Bench
835mm H x 600mm W x 615mm D

Product Code: FPD-08131

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[Shoreline SM150 - 150 Litres Pharmacy Fridge with Solid Door](#)

£798.00

Capacity: 150 Litre
Config: Under Bench
855mm H x 550mm W x 600mm D

Product Code: FPD-04361



[Labcold IntelliCold™ RLDF0510 - 150 Litre Pharmacy Fridge with Solid Door](#)

£860.00

Capacity: 150 Litre
Config: Under Bench
820mm H x 600mm W x 610mm D

Product Code: FPD-09837



[Labcold IntelliCold™ RLDG0510 - 150 Litre Pharmacy Fridge with Glass Door](#)

£1,060.00

Capacity: 150 Litre
Config: Under Bench
820mm H x 600mm W x 610mm D

Product Code: FPD-09838



[LEC PSR151UK - Pharmacy Fridge 151 Litre with Solid Door](#)

£895.00

Capacity: 151 Litre
Config: Under Bench
845mm H x 595mm W x 635mm D

Product Code: FPD-09746



[LEC PGR151UK - Pharmacy Fridge 151 Litre with Glass Door](#)

£1,195.00

Capacity: 151 Litre
Config: Under Bench
845mm H x 595mm W x 635mm D

Product Code: FPD-09755



[Shoreline SM260 - 260 Litres Pharmacy Fridge with Solid Door](#)

£1,098.00

Capacity: 260 Litre
Config: Free Standing
1225mm H x 600mm W x 600mm D

Product Code: FPD-04363



[Shoreline SM260G - 260 Litres Pharmacy Fridge with Glass Door](#)

£1,348.00

Capacity: 260 Litre
Config: Free Standing
1225mm H x 600mm W x 600mm D

Product Code: FPD-04359

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[LEC PSR273UK - Pharmacy Fridge 273
Litre with Solid Door](#)

£1,195.00

Capacity: 273 Litre
Config: Free Standing
1565mm H x 595mm W x 670mm D
Product Code: FPD-09747



[LEC PGR273UK - Pharmacy Fridge 273
Litre with Glass Door](#)

£1,595.00

Capacity: 273 Litre
Config: Free Standing
1565mm H x 595mm W x 660mm D
Product Code: FPD-09756



[Labcold IntelliCold™ RLDG1010 - 300
Litre Pharmacy Fridge with Solid Door](#)

£1,160.00

Capacity: 300 Litre
Config: Free Standing
1500mm H x 600mm W x 660mm D
Product Code: FPD-09839



[Labcold IntelliCold™ RLDG1010 - 300
Litre Pharmacy Fridge with Glass Door](#)

£1,510.00

Capacity: 300 Litre
Config: Free Standing
1500mm H x 600mm W x 600mm D
Product Code: FPD-09840



[LEC PGR334UK - Pharmacy Fridge 334
Litre with Glass Door](#)

£1,995.00

Capacity: 334 Litre
Config: Free Standing
1860mm H x 595mm W x 660mm D
Product Code: FPD-09757



[LEC PSR334UK - Pharmacy Fridge 334
Litre with Solid Door](#)

£1,795.00

Capacity: 334 Litre
Config: Free Standing
1860mm H x 595mm W x 670mm D
Product Code: FPD-09748



[Shoreline SM360G - 360 Litres
Pharmacy Fridge with Glass Door](#)

£1,668.00

Capacity: 360 Litre
Config: Free Standing
1595mm H x 600mm W x 600mm D
Product Code: FPD-04356



[Shoreline SM360 - 360 Litres Pharmacy
Fridge with Solid Door](#)

£1,398.00

Capacity: 360 Litre
Config: Free Standing
1595mm H x 600mm W x 600mm D
Product Code: FPD-04364

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[Labcold IntelliCold™ RLDF1510 - 430 Litre Pharmacy Fridge with Solid Door](#)

£1,960.00

Capacity: 430 Litre
Config: Free Standing
1860mm H x 600mm W x 700mm D

Product Code: FPD-09843



[Labcold IntelliCold™ RLDG1510 - 430 Litre Pharmacy Fridge with Glass Door](#)

£2,125.00

Capacity: 430 Litre
Config: Free Standing
1860mm H x 600mm W x 700mm D

Product Code: FPD-09844



[LEC PE1607C - Pharmacy Fridge 444 Litre with Solid Door](#)

£2,975.00

Capacity: 444 Litre
Config: Free Standing
1688mm H x 685mm W x 910mm D

Product Code: FPD-09749



[LEC PG1607C - Pharmacy Fridge 444 Litre with Glass Door](#)

£3,229.00

Capacity: 444 Litre
Config: Free Standing
1688mm H x 685mm W x 910mm D

Product Code: FPD-09758



[Shoreline SM541 - 544 Litres Pharmacy Fridge with Solid Door](#)

£2,485.00

Capacity: 544 Litre
Config: Free Standing
1745mm H x 750mm W x 700mm D

Product Code: FPD-07642



[Shoreline SM541G - 544 Litres Pharmacy Fridge with Glass Door](#)

£2,850.00

Capacity: 544 Litre
Config: Free Standing
1745mm H x 750mm W x 700mm D

Product Code: FPD-07643



[LEC PSR600UK - Pharmacy Fridge 600 Litre with Solid Door](#)

£3,500.00

Capacity: 600 Litre
Config: Free Standing
1995mm H x 690mm W x 840mm D

Product Code: FPD-10199



[LEC PGR600UK - Pharmacy Fridge 600 Litre with Glass Door](#)

£4,000.00

Capacity: 600 Litre
Config: Free Standing
1995mm H x 690mm W x 840mm D

Product Code: FPD-10201

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[LEC PSR1200UK - Pharmacy Fridge](#)
[1200 Litre with Solid Door](#)

£6,000.00

Capacity: 1200 Litre
Config: Free Standing
1995mm H x 1330mm W x 840mm D
Product Code: FPD-10200



[LEC PGR1200UK - Pharmacy Fridge](#)
[1200 Litre with Glass Door](#)

£6,500.00

Capacity: 1200 Litre
Config: Free Standing
1995mm H x 1330mm W x 840mm D
Product Code: FPD-10202

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

28 Litre Silenco Lockable Pharmacy Absorption Fridge



Product Code: FPD-09818

Dimensions: 580mm H x 422mm W x 393mm D

- The perfect solution for hospitals and health centres to keep heat sensitive medicines between 2-8°C. With magnetic sealing door closing. LED light and acoustic signal as standard. Signal indicates that the door has been opened too long. Signal starts 20 seconds after door opening.
- Internal fan to equalize cabinet temperature and avoid frost formation on the evaporators fin. Trays with transparent front board to keep the inside from warming up by blocking hot air penetration into the cabinet when opening the door. It also displays the min and max temperatures from the last reset.
- Digital thermometer to display the actual inside temperature
- Outside jacket material: polystyrol
- Cooling unit control: thermofuse
- Colour: white
-

Delivery

Carriage

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold IntelliCold™ RLDF0110 - 36 Litre Pharmacy Fridge with Solid Door



Product Code: FPD-09835

Dimensions: 535mm H x 445mm W x 465mm D

- Powered by the unique IntelliCold™, this compact pharmacy and vaccine fridge has been designed from the ground up for the safe storage of cold chain medicines.
- Featuring an all metal interior and exterior, it is built to meet the demands of a professional setting and comes with all the features you would expect on a much larger and expensive professional pharmacy refrigerator.
- Small enough to sit on a bench top or to be wall mounted by way of integral brackets, this versatile pharmacy fridge is ideal for surgery situations, care homes or schools, indeed anywhere where there is a requirement to store limited amount of temperature sensitive vaccines and pharmaceuticals in optimum conditions.
- Configuration: Wall mounted/Bench top
- Operating temp (nominal) 5°C
- Temperature range 2 to 8°C
- Door Type: Solid
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up
- Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Two Shelves



Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PG102C - Pharmacy Fridge 41 Litre with Glass Door



Product Code: FPD-09752

Dimensions: 702mm H x 379mm W x 500mm D

Internal Dimensions: 400mm H x 245mm W x 310mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Small countertop size- For limited space needing medical refrigeration
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Door Finish: Painted aluminium
- Colour: White/green
- Door Type: Glass, Hinged-R/H only, Door Open Depth: 810mm
- Two wire shelves, Number of Shelf Positions: 5, Size: 270x308x6mm
- Automatic defrost
- Climate Class: SN
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Noise Level: <45dB
- Insulation Thickness: 40mm, Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 1.68

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PE102C - Pharmacy Fridge 41 Litre with Solid Door



Product Code: FPD-09742

Dimensions: 702mm H x 379mm W x 500mm D

Internal Dimensions: 400mm H x 245mm W x 310mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Small countertop size - For limited space needing medical refrigeration
- Comes supplied with Wall Bracket
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Door Finish: Painted steel
- Colour: White/green
- Door Type: Solid, Hinged R/H, Door Open Depth: 810Rmm
- Two wire shelves, available positions, Shelf Size: 270x308x6mm
- Automatic defrost
- Climate Class: SN
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Noise Level: <45dB
- Insulation Thickness: 40mm, Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 1.56

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PE109C - Pharmacy Fridge 45 Litre with Solid Door



Product Code: FPD-09743

Dimensions: 530mm H x 470mm W x 465mm D

Internal Dimensions: 400mm H x 410mm W x 270mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- Small countertop size - For limited space needing medical refrigeration
- Comes supplied with Wall Bracket
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Colour: White/green
- Door Type: Solid, Hinged-reversible, Door Open Depth: 910mm
- Two Adjustable Feet at front
- One wire shelf, Size: 402x200x6mm
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Noise Level: <45dB
- Insulation Thickness: 30mm, Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 1.68

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM45 - 45 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-06627

Dimensions: 715mm H x 380mm W x 375mm D

Internal Dimensions: 445mm H x 300mm W x 250mm D

- 45 litres Capacity
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and 2 keys
- CFC & Ammonia free
- Adjustable Shelving x 2
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- 220-230v 50Hz
- Max. Ambient +32°C
- Prices Includes Delivery Mainland UK
- 2 Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP™

53 Litre Silenco Lockable Pharmacy Absorption Fridge



Product Code: FPD-09819

Dimensions: 592mm H x 486mm W x 494mm D

- Digital thermometer to display the actual inside temperature. It also displays the min and max temperatures from the last reset. LED light and acoustic signal as standard. Signal indicates that the door has been opened too long. Signal starts 20 seconds after door opening.
- Internal fan to equalize cabinet temperature and avoid frost formation on the evaporators fin. Trays with transparent front board to keep the inside from warming up by blocking hot air penetration into the cabinet when opening the door.
- Lock/door closing: magnetic sealing
- Cabinet temperature: 2°C to 8°C
- Outside jacket material: polystyrol
- Cooling unit control: thermofuse
- Colour: white

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM60G - 60 Litres Pharmacy Fridge with Glass Door



Product Code: FPD-04358

Dimensions: 660mm H x 470mm W x 455mm D

Internal Dimensions: 525mm H x 375mm W x 310mm D

- Glass Door
- Capacity 60 Litres
- 2 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

Labcold IntelliCold™ RLDF0210 - 66 Litre Pharmacy Fridge with Solid Door



Product Code: FPD-09836

Dimensions: 725mm H x 450mm W x 485mm D

- Powered by the unique IntelliCold™, this compact pharmacy and vaccine fridge has been designed from the ground up for the safe storage of cold chain medicines.
- Featuring an all metal interior and exterior, it is built to meet the demands of a professional setting and comes with all the features you would expect on a much larger and expensive professional pharmacy refrigerator.
- Small enough to sit on a bench top or to be wall mounted by way of integral brackets, this versatile pharmacy fridge is ideal for surgery situations, care homes or schools. Fully lockable and with the added extra of automatic 24/7 integral temperature logging, this pharmacy fridge is designed to be compliant with all current professional guidelines including the DH Green Book, NaTHNaC, GPhC/RPS and RCVS/BSVA.
- Configuration: Bench top/wall mounted
- Door Type: Solid
- Operating temp (nominal): 5°C
- Temperature range: 2 to 8deg;C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Three Shelves



Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM75 - 75 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-07731

Dimensions: 790mm H x 450mm W x 430mm D

Internal Dimensions: 525mm H x 375mm W x 305mm D

- 75 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and 2 keys
- CFC & Ammonia free
- Adjustable Shelving x 2
- Meets RPSGB, NHS, WHO, RCVS performance criteria
- 220-240v 50/60Hz
- Max. Ambient +32°C
- Made in England
- 2 Years Part & Labour Warranty (UK only)
- Price Includes Delivery to Mainland UK

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

LEC PE207C - Pharmacy Fridge 82 Litre with Solid Door



Product Code: FPD-09744

Dimensions: 660mm H x 502mm W x 540mm D

Internal Dimensions: 540mm H x 419mm W x 315mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Door Finish: Painted steel
- Colour: White/green
- Door Type: Solid, Hinged-reversible, Door Open Depth: 976mm
- Two wire shelves, Size: 438x274x6mm
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Noise Level: <43dB
- Insulation Thickness: 30mm
- Energy Consumption kw/24h: 1.46

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PG207C - Pharmacy Fridge 82 Litre with Glass Door



Product Code: FPD-09753

Dimensions: 660mm H x 502mm W x 450mm D

Internal Dimensions: 540mm H x 419mm W x 320mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Door Finish: Painted aluminium
- Colour: White/green
- Door Type: Glass, Hinged-reversible, Door Open Depth: 976mm, Replacement Door Seal
- Open Door Sensor
- Two wire Shelves, Number of Shelf Positions: 2, Size: 438x274x6mm
- Interior Light (Medical)
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Noise Level: <43dB
- Insulation Thickness: 30mm
- Energy Consumption kw/24h: 1.58

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PG307C - Pharmacy Fridge 107 Litre with Glass Door



Product Code: FPD-09754

Dimensions: 830mm H x 495mm W x 620mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted aluminium
- Colour: White/green
- Door Type: Glass, Hinged-reversible, Door Open Depth: 1050mm, Replacement Door Seal
- Adjustable Feet: 1, front
- Two wire shelves, Number of Shelf Positions: 4, Size - 427x374x6mm
- Interior Light (Medical)
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230, Aux Output
- Insulation Thickness: 35mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PE307C - Pharmacy Fridge 107 Litre with Solid Door



Product Code: FPD-09745

Dimensions: 830mm H x 495mm W x 620mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted steel
- Colour: White/green
- Door Type: Solid, Hinged-reversible, Door Open Depth: 1050mm
- Adjustable Feet: 1, front
- Two wire shelves, Number of Shelf Positions: 4, size: 427x374x6mm
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 35mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM161 - 141 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-08131

Dimensions: 835mm H x 600mm W x 615mm D

Internal Dimensions: 667mm H x 441mm W x 439mm D

- Solid Door
- 141 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display with ShoreTouch? controls
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Self-Closing Door with Lock and 2 keys
- High Efficiency - Thick Wall Insulation
- Adjustable Shelving x 3
- Interior Light
- CFC & Ammonia free
- RPSGB, NHS, WHO, RCVS compliant
- 220-230v 50Hz
- Max. Ambient +43°C
- Made in England
- 2 Years Parts & Labour Warranty (UK only)
- Supplied with Certificate of Conformity

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

Shoreline SM161G - 141 Litres Pharmacy Fridge with Glass Door



Product Code: FPD-08130

Dimensions: 835mm H x 600mm W x 615mm D

Internal Dimensions: 667mm H x 441mm W x 439mm D

- Glass Door
- 141 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display with ShoreTouch controls
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Self-Closing Door with Lock and 2 keys
- High Efficiency - Thick Wall Insulation
- Adjustable Shelving x 3
- Interior Light
- CFC & Ammonia free
- RPSGB, NHS, WHO, RCVS compliant
- 220-230v 50Hz
- Max. Ambient +43°C
- Made in England
- 2 Years Parts & Labour Warranty (UK only)
- Supplied with Certificate of Conformity

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

Shoreline SM260 - 260 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-04363

Dimensions: 1225mm H x 600mm W x 600mm D

Internal Dimensions: 1047mm H x 513mm W x 441mm D

- Capacity 260 Litres Gross, 247 Litres Net
- 4 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM150 - 150 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-04361

Dimensions: 855mm H x 550mm W x 600mm D

Internal Dimensions: 665mm H x 460mm W x 385mm D

- Capacity 150 Litres
- 3 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUPTM

Labcold IntelliCold™ RLDG0510 - 150 Litre Pharmacy Fridge with Glass Door



Product Code: FPD-09838

Dimensions: 820mm H x 600mm W x 610mm D



- A professional 150 litre under bench pharmacy refrigerator powered by the Labcold IntelliCold™ controller. Designed for the safe storage of cold chain pharmaceuticals, this refrigerator assists busy professionals meet the regulatory requirements with an easy to operate intuitive four button controller and 13 month temperature data storage which can be downloaded to any PC or Mac. For ease of viewing the contents and to minimize door openings, the light can be set to operate independently of the door and even left on all the time.
- Reliable and robust, this refrigerator is engineered to meet the demands of a busy practice, surgery or pharmacy and is backed by the Labcold FREE 5 year warranty.
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Underbench
- Operating temp (nominal): 5°C
- Temperature range: 2 to 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Three Shelves
- Door Type: Glass
- Access port

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PSR273UK - Pharmacy Fridge 273 Litre with Solid Door



Product Code: FPD-09747

Dimensions: 1565mm H x 595mm W x 670mm D

Internal Dimensions: 1382mm H x 480mm W x 445mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Anti-microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Interior Xw light - For better visibility of contents
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Type: Solid, Hinged-reversible, Door Open Depth: 1185mm
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Five wire Shelves, Number of Shelf Positions: 12
- Interior Light (Medical)
- Automatic defrost
- Climate Class: SN
- Refrigerant Type: R600a
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 220-240v/50hz, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold IntelliCold™ RLDF0510 - 150 Litre Pharmacy Fridge with Solid Door



Product Code: FPD-09837

Dimensions: 820mm H x 600mm W x 610mm D



- A professional 150 litre under bench pharmacy refrigerator powered by the Labcold IntelliCold™ controller. Designed for the safe storage of cold chain pharmaceuticals, this refrigerator helps busy professionals meet regulatory requirements laid down in the DH Green Book with an easy to operate intuitive four button controller and 13 month temperature data storage. This data can be downloaded to a free SD card (available when the warranty is registered) and stored on a PC or Mac. Reliable and robust, this refrigerator is engineered to meet the demands of a busy practice or pharmacy and is backed by the Labcold FREE 5 year parts and labour warranty.
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Underbench
- Door Type: Solid
- Operating temp (nominal): 5°C
- Temperature range: 2 - 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Three Shelves
- Access port

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PGR151UK - Pharmacy Fridge 151 Litre with Glass Door



Product Code: FPD-09755

Dimensions: 845mm H x 595mm W x 635mm D

Internal Dimensions: 705mm H x 500mm W x 450mm D

- PSR151 solid door refrigerator is ideal for the under-counter storage of vaccines and comes complete with our new advanced controller and a fitted lock for added security
- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Antimicrobial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening.
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted steel
- Door Type: Solid, Hinged-reversible, Door Open Depth: 1165mm
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Three wire shelves, Number of Shelf Positions: 6, Size: 400x520x6mm
- Automatic defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 45mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PGR273UK - Pharmacy Fridge 273 Litre with Glass Door



Product Code: FPD-09756

Dimensions: 1565mm H x 595mm W x 660mm D

Internal Dimensions: 1382mm H x 480mm W x 445mm D

- PGR273 glass door refrigerator is ideal for the storage of vaccines and comes complete with our new advanced controller and a fitted lock for added security.
- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Interior Xw light - For better visibility of contents without having to open the door
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted aluminium
- Door Type: Glass, Hinged-reversible, Door Open Depth: 1180mm, Replacement Door Seal
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Five wire Shelves, Number of Shelf Positions: 12
- Automatic defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 220-240v/50hz, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PSR151UK - Pharmacy Fridge 151 Litre with Solid Door



Product Code: FPD-09746

Dimensions: 845mm H x 595mm W x 635mm D

Internal Dimensions: 705mm H x 500mm W x 450mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Antimicrobial handle
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening.
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted steel
- Colour: White/green
- Door Type: Solid, Hinged-reversible, Door Open Depth: 1165mm
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Three wire Shelves, Number of Shelf Positions: 6, size: 400x520x6mm
- Automatic defrost
- Climate Class: SN
- Refrigerant Type: R600a
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 45mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM260G - 260 Litres Pharmacy Fridge with Glass Door



Product Code: FPD-04359

Dimensions: 1225mm H x 600mm W x 600mm D

Internal Dimensions: 1047mm H x 513mm W x 441mm D

- Glass Door
- Interior Light
- Capacity 260 Litres Gross, 252 Litres Net
- 4 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP_™

Labcold IntelliCold™ RLDG1010 - 300 Litre Pharmacy Fridge with Glass Door



Product Code: FPD-09840

Dimensions: 1500mm H x 600mm W x 600mm D

- With a toughened safety glass door, this 300 litre pharmacy fridge, fitted with the Labcold IntelliCold™ controller, is purposed designed for the safe storage of temperature sensitive vaccines and medicines. With an external light switch that can operate independently of the fridge door, this model can not only form an attractive display feature but allows for quick visual stock takes without the need for opening the door. Backed by the Labcold FREE 5 year parts and labour warranty, there has never been a better way to store cold chain pharmaceuticals.
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Upright
- Door Type: Glass
- Operating temp (nominal): 5°C
- Temperature range: 2 - 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Six Shelves
- Access port

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.



FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold IntelliCold™ RLDF1010 - 300 Litre Pharmacy Fridge with Solid Door



Product Code: FPD-09839

Dimensions: 1500mm H x 600mm W x 660mm D

- A freestanding 300 litre pharmacy refrigerator fitted with the Labcold IntelliCold™ controller to ensure the safe storage of cold chain pharmaceuticals. The easy to operate push button controller clearly displays the internal temperature and displays minimum and maximum temperature since the last temperature re-set at the touch of the button. Immediate high/low temperature alarms and controller battery back up which means it functions even in the event of a power failure, makes this fridge one of the safest for vaccine storage. For extra security, this refrigerator is backed by the Labcold FREE 5 year parts and labour warranty.
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Upright
- Door Type: Solid
- Operating temp (nominal): 5°C
- Temperature range: 2 to 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Six Shelves
- Access port



Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PG1607C - Pharmacy Fridge 444 Litre with Glass Door



Product Code: FPD-09758

Dimensions: 1688mm H x 685mm W x 910mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Interior 12w light - For better visibility of contents without having to open the door
- Front lockable roller castors - To make the fridge easier to move for cleaning and maintenance, and secure it in place
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted aluminium
- Colour: White/green
- Door Type: Glass, Hinged-R/H only, Replacement Door Seal
- Five wire shelves, Number of Shelf Positions: 30, Size: 536x609x6mm
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 40mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PGR334UK - Pharmacy Fridge 334 Litre with Glass Door



Product Code: FPD-09757

Dimensions: 1860mm H x 595mm W x 660mm D

Internal Dimensions: 1682mm H x 480mm W x 445mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Interior Xw light - For better visibility of contents without having to open the door
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted aluminium
- Colour: White/green
- Door Type: Glass, Hinged-reversible, Door Open Depth: 1180mm, Replacement Door Seal
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Six wire shelves, Number of Shelf Positions: 14, Size: 385x500x6mm
- Automatic defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 220-240v/50hz, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PSR334UK - Pharmacy Fridge 334 Litre with Solid Door



Product Code: FPD-09748

Dimensions: 1860mm H x 595mm W x 670mm D

Internal Dimensions: 1682mm H x 480mm W x 445mm D

- Our large capacity PSR334 solid door refrigerator is ideal for the storage of vaccines and comes complete with our new advanced controller and a fitted lock for added security.
- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Antimicrobial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening.
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Interior Xw light - For better visibility of contents
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Type: Solid, Hinged-reversible, Door Open Depth: 1185mm
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Six wire Shelves, Number of Shelf Positions: 14, Size: 6: 385x500x6mm
- Interior Light (Medical)
- Automatic defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 220-240v/50hz, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PE1607C - Pharmacy Fridge 444 Litre with Solid Door



Product Code: FPD-09749

Dimensions: 1688mm H x 685mm W x 910mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Front lockable roller castors - To make the fridge easier to move for cleaning and maintenance, and secure it in place
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Type: Solid, Hinged-R/H only, Replacement Door Seal
- Five wire shelves, Number of Shelf Positions: 30, Size: 536x609x6mm
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 40mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM360G - 360 Litres Pharmacy Fridge with Glass Door



Product Code: FPD-04356

Dimensions: 1595mm H x 600mm W x 600mm D

Internal Dimensions: 1418mm H x 513mm W x 441mm D

- Glass Door
- Interior Light
- Gross Volume: 352 Litres, Net Volume 345 Litres
- 5 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48
HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP_™

Shoreline SM360 - 360 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-04364

Dimensions: 1595mm H x 600mm W x 600mm D

Internal Dimensions: 1418mm H x 513mm W x 441mm D

- Capacity 360 Litres Gross, 352 Litres Net
- 5 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM541G - 544 Litres Pharmacy Fridge with Glass Door



Product Code: FPD-07643

Dimensions: 1745mm H x 750mm W x 700mm D

Internal Dimensions: 1452mm H x 600mm W x 560mm D

- High Efficiency Model
- 544 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Remote Alarm Terminals
- Automatic Defrosting
- Door Lock and 2 keys
- Castors (2 front locking)
- Self-Closing Door
- High Efficiency Model
- Thick Wall Insulation
- 1.4 kWhr/24hr
- Adjustable Shelving x 4
- Interior Light
- CFC & Ammonia free
- RPSGB, NHS, WHO, RCVS compliant
- 220-230v 50Hz
- Max. Ambient +38°C
- Made in England
- 2 Years Part & Labour Warranty (UK only)
- Price Includes Delivery to Mainland UK

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

Labcold IntelliCold™ RLDG1510 - 430 Litre Pharmacy Fridge with Glass Door



Product Code: FPD-09844

Dimensions: 1860mm H x 600mm W x 700mm D

- A 430 litre glass door pharmacy refrigerator powered by the Labcold IntelliCold™ controller. Two internal fans are fitted to this refrigerator, giving it superior temperature uniformity making it one of the safest ways to store large quantities of cold chain pharmaceuticals. The large safety glass door reduces the need for door openings as the internal light can be operated from the exterior of the fridge so can be turned on independently of the door opening. The easy four button operation, plus the capacity to store temperature data for up to 13 months assist with regulatory compliance and for extra peace of mind it is backed by the Labcold FREE 5 year parts and labour warranty.
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Upright
- Door Type: Glass
- Operating temp (nominal): 5°C
- Temperature range: 2 to 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- FREE SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Six Shelves
- Access port



Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold IntelliCold™ RLDF1510 - 430 Litre Pharmacy Fridge with Solid Door



Product Code: FPD-09843

Dimensions: 1860mm H x 600mm W x 700mm D

- A slim line 430 litre pharmacy refrigerator powered by the Labcold IntelliCold™ controller and fitted with two fans for the optimum in temperature stability. Designed exclusively for the safe storage of vaccines and other temperature sensitive pharmaceuticals, this fridge makes compliance with regulatory guidelines simple by featuring an intuitive four button controller and 13 month temperature data logging. A free SD card for temperature data downloads is available when you register for your FREE 5 year parts and labour warranty!
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Upright
- Door Type: Solid
- Operating temp (nominal): 5°C
- Temperature range: 2 to 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Six Shelves
- Access port



Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited
The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM541 - 544 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-07642

Dimensions: 1745mm H x 750mm W x 700mm D

Internal Dimensions: 1452mm H x 600mm W x 560mm D

- High Efficiency Model
- 544 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Remote Alarm Terminals
- Automatic Defrosting
- Door Lock and 2 keys
- Castors (2 front locking)
- Self-Closing Door
- High Efficiency Model
- Thick Wall Insulation
- 1.0 kWhr/24hr
- Adjustable Shelving x 5
- CFC & Ammonia free
- RPSGB, NHS, WHO, RCVS compliant
- 220-230v 50Hz
- Max. Ambient +43°C
- Made in England
- 2 Years Part & Labour Warranty (UK only)
- Price Includes Delivery to Mainland UK

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

LEC PSR600UK - Pharmacy Fridge 600 Litre with Solid Door



Product Code: FPD-10199

Dimensions: 1995mm H x 690mm W x 840mm D

Internal Dimensions: 1386mm H x 560mm W x 653mm D

- Conforms to MHRA guidelines.
- NEW Advanced Controller
- Records temperatures accurate to 0.1°C
- Fitted lock as standard maintains security and ensures the unit can only be accessed by approved personnel
- Fitted Castors for ease of movement
- Wire shelves and basket enables optimum air flow inside the unit
- Temperature range +2°C to +8°C
- Colour White/green
- External finish White
- Internal finish Stainless steel
- Type of plastic SRL800 (HIPS/ABS Mix)
- Door finish Painted stainless steel
- Door type Solid, Hinged
- Door open depth 1145mm
- Adjustable feet (number & position) Fitted castors
- Temperature display on panel
- Temperature alarm
- Five shelves, 22 available positions
- Shelf size 645mm x 523mm
- Shelf material Plastic coated wire
- One plastic coated wire baskets
- Fan
- Automatic defrost
- External digital controller
- Climate class N
- Refrigerant type R134a
- CFC, HFC, HCFC & Ammonia free
- UK 3 pin plug, Other plug types available, Voltage supply 230V, Aux output
- Insulation thickness: 60mm, Insulation type: Cyclopentane
- Optional drawer system available

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PGR600UK - Pharmacy Fridge 600 Litre with Glass Door



Product Code: FPD-10201

Dimensions: 1995mm H x 690mm W x 840mm D

Internal Dimensions: 1386mm H x 560mm W x 653mm D

- Conforms to MHRA guidelines
- NEW Advanced Controller
- Records temperatures accurate to 0.1°C
- Fitted lock as standard maintains security and ensures the unit can only be accessed by approved personnel
- Fitted castors for ease of movement
- Wire Shelves and Basket enables optimum air flow inside the unit
- Temperature range +2°C to +8°C
- Colour White/green
- External finish White
- Internal finish Stainless steel
- Type of plastic SRL800 (HIPS/ABS Mix)
- Door finish Painted stainless steel
- Door type Glass, Hinged
- Door open depth 1145mm
- Adjustable feet (number & position) Fitted castors
- Temperature display on panel
- Temperature alarm
- Five shelves with 22 available positions
- Shelf size 645mm x 523mm
- One plastic coated wire basket
- Interior light
- Fan
- Automatic defrost
- External digital controller
- Climate class N
- Refrigerant type R134a
- CFC, HFC, HCFC & Ammonia free
- UK 3 pin plug, Other plug types available, Voltage supply 230V, Aux output
- Insulation thickness 60mm, Insulation type Cyclopentane
- Optional drawer system available

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PGR1200UK - Pharmacy Fridge 1200 Litre with Glass Door



Product Code: FPD-10202

Dimensions: 1995mm H x 1330mm W x 840mm D

Internal Dimensions: 1386mm H x 1220mm W x 653mm D

- Conforms to MHRA guidelines
- NEW Advanced Controller records temperatures accurate to 0.1°C
- Fitted Lock as Standard maintains security and ensures the unit can only be accessed by approved personnel
- Fitted Castors for ease of movement
- Wire Shelves and Basket enables optimum air flow inside the unit
- Temperature range +2°C to +8°C
- Colour White/green
- External finish White
- Internal finish Stainless steel
- Type of plastic SRL800 (HIPS/ABS Mix)
- Door finish Painted stainless steel
- Door type Glass, Hinged
- Door open depth 1145mm
- Adjustable feet (number & position) Fitted castors
- Temperature display on panel
- Temperature alarm
- Ten shelves 44 positions available
- Shelf size 645mm x 523mm
- Two plastic coated wire baskets
- Interior light
- Fan
- Automatic defrost
- External digital controller
- Climate class N
- Refrigerant type R134a
- CFC free, HFC, HCFC & Ammonia free
- UK 3 pin plug, Other plug types available, Voltage supply 230V, Aux output
- Insulation thickness 60mm, Insulation type Cyclopentane
- Optional drawer system available

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PSR1200UK - Pharmacy Fridge 1200 Litre with Solid Door



Product Code: FPD-10200

Dimensions: 1995mm H x 1330mm W x 840mm D

Internal Dimensions: 1386mm H x 1220mm W x 653mm D

- Conforms to MHRA guidelines
- NEW Advanced Controller
- Records temperatures accurate to 0.1°C
- Fitted lock as standard maintains security and ensures the unit can only be accessed by approved personnel
- Fitted castors for ease of movement
- Wire shelves and basket enables optimum air flow inside the unit
- Temperature range +2°C to +8°C
- Colour White/green
- External finish White
- Internal finish Stainless steel
- Type of plastic SRL800 (HIPS/ABS Mix)
- Door finish Painted stainless steel
- Door type Solid, Hinged
- Door open depth 1145mm
- Adjustable feet (number & position) Fitted castors
- Temperature display on panel
- Temperature alarm
- Ten shelves, 44 positions available
- Shelf size 645mm x 523mm
- Two plastic coated wire baskets
- Fan
- Automatic defrost
- External digital controller
- Climate class N
- Refrigerant type R134a
- CFC, HFC, HCFC & Ammonia
- UK 3 pin plug, Other plug types available, Voltage supply 230V, Aux output
- Insulation thickness 60mm, Insulation type Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Advanced Pharmacy Fridges

The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.



[Labcold RPF05043 - 150 Litre
Advanced Pharmacy Fridge with Solid
Door](#)

£1,950.00

Capacity: 150 Litre
Config: Under Bench
820mm H x 600mm W x 620mm D

Product Code: FPD-07551



[Labcold RPF05043 - 150 Litre
Advanced Pharmacy Fridge with Glass
Door](#)

£2,145.00

Capacity: 150 Litre
Config: Under Bench
820mm H x 600mm W x 620mm D

Product Code: FPD-07552



[Labcold RPF21042 - 571 Litre
Advanced Pharmacy Fridge with Solid
Door](#)

£3,555.00

Capacity: 571 Litre
Config: Free Standing
1980mm H x 720mm W x 820mm D

Product Code: FPD-07553



[Labcold RPF21042 - 571 Litre
Advanced Pharmacy Fridge with Glass
Door](#)

£4,100.00

Capacity: 571 Litre
Config: Free Standing
1980mm H x 720mm W x 820mm D

Product Code: FPD-07554



[Labcold RPF44042 - 1245 Litre
Advanced Pharmacy Fridge with Solid
Door](#)

£5,425.00

Capacity: 1245 Litre
Config: Free Standing
1980mm H x 1440mm W x 820mm D

Product Code: FPD-06566



[Labcold RPF44042 - 1245 Litre
Advanced Pharmacy Fridge with Glass
Door](#)

£6,305.00

Capacity: 1245 Litre
Config: Free Standing
1980mm H x 1440mm W x 820mm D

Product Code: FPD-07555

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPF05043 - 150 Litre Advanced Pharmacy Fridge with Solid Door



Product Code: FPD-07551

Dimensions: 820mm H x 600mm W x 620mm D

Internal Dimensions: 680mm H x 540mm W x 350mm D

- Robustly constructed from Stainless Steel this refrigerator from the Labcold Advanced range is designed to provide close temperature control for temperature sensitive products.
- The fan circulated interior has no internal 'step', making it easy to clean and stock and to allow for maximum air flow for even temperature distribution. The two shelves are fully adjustable, so it is possible to maximise the space available.
- Insulated with Envirofoam, which is derived from plant rather than petrochemical oils, the refrigerant is also ammonia, HCFC and CFC free.
- Fitted with a sturdy door lock, this Advanced refrigerator is ideal for busy environments such as busy laboratories or pharmacies.
- Config: Under Counter
- Energy Consumption 3.75 kWh/24h
- Operating Temperature +4°C
- Temperature Set Point Range +2°C to +8°C
- Interior & Exterior Construction Stainless Steel
- Two Shelves
- Shelf dimensions (W xD) 530 x 325mm
- Door Hinging Right hand only
- Door Style: Solid
- Internal Light Glass Door Model only
- Rear Castors
- Door Lock
- Access Port Optional Extra
- Digital Temperature Control, Digital Temperature Display
- High Temperature alarm, Low Temperature alarm, Door Open alarm
- Auto defrost
- Fan Assisted Cooling
- Environmentally friendly insulation
- Eco Friendly Refrigerant
- Ammonia Free
- DH Green Book, RPS/GPhC, NaTHNaC and RCVS Guidelines compliant
- Voltage Supply: 230V/50Hz, UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPF21042 - 571 Litre Advanced Pharmacy Fridge with Glass Door



Product Code: FPD-07554

Dimensions: 1980mm H x 720mm W x 820mm D

- The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.
- Configuration: Free-standing
- Capacity: 571 Litre
- Auto Defrost
- Shelves: 4
- Operating temp (nominal) +4°C
- Door Type: Glass
- Lockable castors
- Precise microprocessor temperature control
- Digital temperature display and alarm
- Automatic defrost
- Adjustable shelving
- Lockable castors
- Lockable door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPF21042 - 571 Litre Advanced Pharmacy Fridge with Solid Door



Product Code: FPD-07553

Dimensions: 1980mm H x 720mm W x 820mm D

- The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.
- Configuration: Free-standing
- Capacity: 571 Litre
- Auto Defrost
- Shelves: 4
- Operating temp (nominal) +4°C
- Door Type: Solid
- Lockable castors
- Precise microprocessor temperature control
- Digital temperature display and alarm
- Automatic defrost
- Adjustable shelving
- Lockable castors
- Lockable door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPF21042 - 571 Litre Advanced Pharmacy Fridge with Glass Door



Product Code: FPD-07554

Dimensions: 1980mm H x 720mm W x 820mm D

- The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.
- Configuration: Free-standing
- Capacity: 571 Litre
- Auto Defrost
- Shelves: 4
- Operating temp (nominal) +4°C
- Door Type: Glass
- Lockable castors
- Precise microprocessor temperature control
- Digital temperature display and alarm
- Automatic defrost
- Adjustable shelving
- Lockable castors
- Lockable door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPFR44042 - 1245 Litre Advanced Pharmacy Fridge with Solid Door



Product Code: FPD-06566

Dimensions: 1980mm H x 1440mm W x 820mm D

- The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.
- Configuration: Free-standing
- Capacity: 1245 Litre
- Shelves: 8
- Operating temp (nominal): +4°C
- Door Type: Solid
- Lockable castors
- Precise microprocessor temperature control
- Digital temperature display and alarm
- Automatic defrost
- Adjustable shelving
- Lockable castors
- Lockable door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPF44042 - 1245 Litre Advanced Pharmacy Fridge with Glass Door



Product Code: FPD-07555

Dimensions: 1980mm H x 1440mm W x 820mm D

- The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.
- Configuration: Free-standing
- Capacity: 1245 Litre
- Shelves: 8
- Operating temp (nominal): +4°C
- Door Type: Double Glass
- Lockable castors
- Precise microprocessor temperature control
- Digital temperature display and alarm
- Automatic defrost
- Adjustable shelving
- Lockable castors
- Lockable door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Delivery

48 Hour delivery is available on most refrigerators through-out the UK mainland, larger units may require a longer lead time, please contact our sales staff whilst ordering.

Approximate Delivery lead times can be confirmed during our quotation process, typical delivery time is 10-14 days from receipt of payment.

**Drivers will not carry goods into buildings.

Goods normally arrive by truck on a pallet and are delivered Kerbside. Please note you may need lifting equipment or a number of strong workers to handle the order when it arrives. If you have any special requirements please [contact us](#).

Report any damages or shortages immediately, please read delivery notes carefully and use contact numbers provided or call FPD Group Ltd.

For more detail refer to [FPD Group Ltd Terms and Conditions](#) documentation.

**Large Products will require a two person delivery and installation, this is built into the price, the larger fridges also may require a longer delivery time of 3-4 weeks please check whilst ordering.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP[™]

How To Order

NHS Purchase Order, MOD Purchase Order, Local Authority Purchase Order

All official government, NHS, MOD, Police, Educational and Local Authority purchase orders are accepted. Fax your purchase order to 0115 727 0123 or email to info@pharmacy-fridge.co.uk.

Bank Transfer

Sent bank transfers to:

Payee	FPD Group Limited		
Bank	HSBC Bank plc		
Sortcode	401903	IBAN	GB97MIDL40190302074613
Account Number	02074613	Swift	MIDLGB2104A

If possible send confirmation of your transaction, to speed up the order process.

Cheques

Make cheques payable to FPD Group Limited and post to:

FPD Group Limited
The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Credit/Debit Card

Payment from UK based businesses can be taken over the telephone or processed via the website [shopping cart](#). Customers with issued quotations may process their order online using our [Purchase Code](#) web page.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Pharmacy Fridge Exports

Since 1961 FPD Group Ltd's family have been supplying medical equipment. Our extensive experience and technical knowledge now concentrates on providing clients with a quality selection of products for hospitals and healthcare organisations worldwide. We receive export enquiries and orders from all parts of the world and have established customers in UK, Europe, Africa and Asia. Export business activities are supervised by either our UK or Asia offices.

Our client base ranges from overseas ministry of health, defence, hospitals, surgeries, clinics and doctors to aid agencies and disaster relief. FPD Group Ltd export almost any medical product listed on our websites.

FPD estimators have the expertise to supply a single item or to fully equip many types of healthcare facilities. As well as the initial supply of equipment we also offer room design services using our in-house Autodesk Revit Architecture 2013 CAD systems.

** Export items are normally subject to minimum order values, additional protective packing charges and export shipping costs. Prices shown are Ex Works unless otherwise agreed in writing.*



FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

About Pharmacy Fridge

Pharmacy Fridge supply a wide range of fridges meeting the specifications of the healthcare industry. To buy simply add your required product to the shopping cart and submit your order.

For dispensary, surgery, ward and general use for the storage of drugs and vaccines under controlled conditions in line with the "Code of Ethics and Standards" of the Royal Pharmaceutical Society of Great Britain. Where necessary to achieve the required temperature control internal air-circulating fans are fitted.

Temperature display, alarms, max/min temperature recording, and door locks are standard and all models are CFC, HCFC and Ammonia free.

Are you storing your vaccines correctly? Recent research has shown an alarming level of use of domestic refrigerators to store vaccines. Did you know that not only can this damage the effectiveness of the vaccines it also poses a security risk for medicines that should be locked away?

FPD Group Ltd

FPD Group Limited is a family business with over 50 years experience in the supply and installation of a wide range of Healthcare Equipment, Medical Refrigeration and Hygienic Storage Furniture. Our customer base includes organisations such as NHS Trusts, PCT's, HMP Prison service, Ministry of Defence, GP Practices, Pharmacy Shops and Dispensaries, Universities, Schools, Dental Practices, Laboratories, Veterinary Practices, Care Homes and Private Buyers. We also export to the Irish Republic, Asia, United Arab Emirates and Qatar. In the UK the Company is based in Exeter, Devon with an additional office in Asia.

Our product range includes pharmacy storage solutions, refrigeration, controlled drugs storage, medical trolleys, clinical waste bins, medical lighting solutions, medical and healthcare furniture ranges and a number of bespoke solutions on request.

FPD Group has an in-house CAD visualisation department and accepts most file formats from your architects' plans, please ask for a quotation.

We are fully committed to the provision of a comprehensive range of professional products, services and furniture. Our philosophy is to provide a complete and appropriate solution to our customers' requirements.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Contact Pharmacy Fridge

FPD Group Ltd have over 40 years experience in the healthcare and shopfitting industry. We offer consultancy, design and installation services. Along with our Pharmacy Shelving and Storage solutions we also provide a wide range of medical equipment and healthcare furniture. Please visit www.fpd-group.co.uk for more information on our medical products.

For general enquiries complete the form below and we will contact you ASAP, please provide as much information as possible including dimensions where required. Please email any plans or photos of installations directly to sales@medicine-fridge.co.uk.

Name	<input type="text"/>	Address	<input type="text"/>
Company	<input type="text"/>	...	<input type="text"/>
Telephone	<input type="text"/>	City	<input type="text"/>
Fax	<input type="text"/>	State	<input type="text"/>
Email	<input type="text"/>	PostalCode	<input type="text"/>
		Country	<input type="text"/>
Comments	<input type="text"/>		

Submit Enquiry

Sales

Tel: 0113 350 8696
Fax: 0115 727 0123
Email: sales@pharmacy-shelving.co.uk

Payments & Billing

Tel: 01392 927 567

Address

FPD Group Limited
The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk